

Initial S™: Single use sterile kits for hand surgery

Initial S™ - Hand - Instrumentation kit content*

STERILE INSTRUMENTATION FOR Ø2.25 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC2.25H	Kit for Ø2.25 mm cannulated screws - Hand	Ø1.7 mm quick coupling drill bit - cannula Ø0.9 mm - L93mm	1
		3 in 1 instrument for Ø2.25 mm screws	1
		Double guide Ø0.9 mm / Ø1.8 mm for Ø2.25 mm screws	1
		5.8 mm single use handle	1
		Pin Ø0.8 L150 mm	2
		Pin Ø0.8 L120 mm	1

STERILE INSTRUMENTATION FOR Ø2.6 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC2.6H	Kit for Ø2.6 mm cannulated screws - Hand	Ø2.0 mm quick coupling drill bit - cannula Ø1.1 mm - L93 mm	1
		3 in 1 instrument for Ø2.6 mm screws	1
		Double guide Ø1.1 mm / Ø2.1 mm for Ø2.6 mm screws	1
		5.8 mm single use handle	1
		Pin Ø1.0 L150 mm	2
		Pin Ø1.0 L120 mm	1

STERILE INSTRUMENTATION FOR Ø3.0 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC3.0H	Kit for Ø3.0 mm cannulated screws - Hand	Ø2.2 mm quick coupling drill bit - cannula Ø1.1 mm - L93 mm	1
		3 in 1 instrument for Ø3.0 mm screws	1
		Double guide Ø1.1 mm / Ø2.3 mm for Ø3.0 mm screws	1
		5.8 mm single use handle	1
		Pin Ø1.0 L150 mm	2
		Pin Ø1.0 L120 mm	1


Indications :

A kit of the INITIAL S™ range consists of a set of instruments or a set of instruments and implants (washers), to be used only in trauma or orthopaedic surgery, when NewclipTechnics screws for osteosynthesis are fitted. Implants of the Stand-Alone Screws range are intended for fractures fixation, osteotomies and arthrodeses of bones in adults and adolescents, appropriate for the size of the device.

Contra-indications :

- Serious vascular deterioration, bone devitalization,
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implants into the bone,
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

*Sterile instrumentation for use with Stand-Alone Screws - Handmotion range. Available separately.



Example of instrumentation kit content for Ø3.0 mm screws.



NEWCLIP-TECHNICS




Initial S™: Single use sterile kits for hand surgery

Initial S™ - Hand - Implant references*


Ø2.25 MM CANNULATED SCREWS

Ref.	Description
H0.9IFT2.25L14-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 14 mm - STERILE
H0.9IFT2.25L16-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 16 mm - STERILE
H0.9IFT2.25L18-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 18 mm - STERILE
H0.9IFT2.25L20-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 20 mm - STERILE
H0.9IFT2.25L22-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 22 mm - STERILE
H0.9IFT2.25L24-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 24 mm - STERILE
H0.9IFT2.25L26-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 26 mm - STERILE
H0.9IFT2.25L28-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 28 mm - STERILE
H0.9IFT2.25L30-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 30 mm - STERILE
H0.9IFT2.25L32-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 32 mm - STERILE
H0.9IFT2.25L34-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 34 mm - STERILE




Ø2.6 MM CANNULATED SCREWS

Ref.	Description
H1.1IFT2.6L14-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 14 mm - STERILE
H1.1IFT2.6L16-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 16 mm - STERILE
H1.1IFT2.6L18-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 18 mm - STERILE
H1.1IFT2.6L20-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 20 mm - STERILE
H1.1IFT2.6L22-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 22 mm - STERILE
H1.1IFT2.6L24-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 24 mm - STERILE
H1.1IFT2.6L26-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 26 mm - STERILE
H1.1IFT2.6L28-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 28 mm - STERILE
H1.1IFT2.6L30-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 30 mm - STERILE
H1.1IFT2.6L32-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 32 mm - STERILE
H1.1IFT2.6L34-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 34 mm - STERILE



Ø3.0 MM CANNULATED SCREWS

Ref.	Description
H1.1IFT3.0L14-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 14 mm - STERILE
H1.1IFT3.0L16-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 16 mm - STERILE
H1.1IFT3.0L18-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 18 mm - STERILE
H1.1IFT3.0L20-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 20 mm - STERILE
H1.1IFT3.0L22-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 22 mm - STERILE
H1.1IFT3.0L24-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 24 mm - STERILE
H1.1IFT3.0L26-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 26 mm - STERILE
H1.1IFT3.0L28-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 28 mm - STERILE
H1.1IFT3.0L30-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 30 mm - STERILE
H1.1IFT3.0L32-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 32 mm - STERILE
H1.1IFT3.0L34-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 34 mm - STERILE



NEWCLIP TECHNICS
PA de la Lande Saint Martin - 45 rue des Garotières
44115 Haute Goulaine (France)
P: +33 (0)2 28 21 23 25 - Fax: +33 (0)2 40 63 68 37
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP USA
642 Larkfield Center
Santa Rosa CA 95403, USA
P: + 1 707 230 5078
customerservice@newclipusa.com
www.newclipusa.com

NEWCLIP GmbH
Pröllstraße 11,
D-86157 Augsburg, Deutschland
P: +49 (0)821 650 749 40
info@newclipgmbh.com
www.newclipgmbh.de

NEWCLIP Australia
3B/11 Donkin Street
West E,d 4101, Australie
P: +61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP TECHNICS Japan K.K
KKK Bldg 502, 3-18-1 Asakusabashi,
Taito-Ku, Tokyo, 111-0053 Japan
P: +81 (0)3 58 25 49 81
Fax: +81 (0)3 58 25 49 86
www.newcliptechnics.fr

The information presented in this brochure is intended to demonstrate a Newclip Technics product. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.